


Scan

Automatic colony counters & inhibition zone readers



interscience



interscience

Our quality for your lab

- Designer and manufacturer for microbiological analyses
- Made in France
- Colony counting specialist
- Worldwide distribution network

Scan

High quality analyses, full traceability

Scan 300, Scan 500 & Scan 1200

High technology automatic colony counters


With a digital camera and high technology software, they can be linked to a PC via a USB connection. They count all colonies on a Petri dish in less than 1 second and provide a complete, fast, accurate and traceable reading of the result.

■ Bacterial enumeration

- Food analyses
- Total flora analyses
- Aerobic & anaerobic bacterial enumeration, yeasts, lactobacillus...
- Pathogenic bacteria research
- Environmental research
- Pharmaceutical analyses
- Medical analyses
- Cosmetics analyses

■ Inhibition zones

- **Pharmaceutical industry, medical research & hospitals**
(antibiograms, resistance tests to pathogenic microbes, medical diagnoses...)
- **Food industry**
(Tests on lactic ferments & for dairy ingredients industry...)




Scan

- Automatic colony counters: Easy settings
- Inhibition zone readers*
- Data traceability and full report









High performance

- > Count colonies of numerous media
- > Reading of chromogenic media⁽¹⁾: Colored differentiation of colonies (up to 4 different colors on the same dish)
- > Inhibition zone measurement



Live image

- > Fits any type of dish: automatic adjustment of contrast and lighting
- > High-definition color image
- > Each colony is marked with a cross
- > Powerful zoom: up to x28

Display details			
Sample N°	Count	Dilution	
 E COLI	174	1 2.0	
 E COLI	353	1 5.0	
 SPIRAL	47	1/1000 9.0	
 PETRIFILM	89	1 8.0	
 RIDA COUNT	179	1 1.0	
 Filtering Membrane	111	1 1.0	

Instant results

- > Up to 1000 colonies detected in 1 second
- > Counts 30 dishes in 5 minutes (in real condition with presetting)
- > Reproducible and standardized results
- > **Scan** results: instant and automatic



(1) on Scan 500 & Scan 1200





Easy-to-use

- Counting in 1 click <
- All functions in 1 single window <
- Custom parameters: day, users, project... <




Dark Field technology

- Display of every colony <
- Optimized lighting & contrast <
- Long lasting LED lighting <
- 6 lighting combinations <



Traceability & reporting

- Automatic archiving and printing of data: < pictures, comments & results
- Export to EXCEL™, PDF, JPEG, PNG, BMP <
- Barcode reader <
- Connection to LIMS network <



www.interscience.com




(2) Free update of the software during warranty period (3) After registration of the warranty form

Efficiency & time saving


Instant results

Thanks to the live image display of the Petri dish on your computer, count more than **1000 CFU/s** on all media. Each counted colony is marked with a cross and the result is automatically saved.



Easy settings

Choose your pre-set parameters for Petri dishes:




Also available on **Scan 1200**:

- MC-Media Pads™: AC, CC, EC/CC, SA
- Petrifilm™: AC, ETB, CC, EC/CC, EC
- Compact Dry™: TC, CF, EC, ETB
- Easy Plate™
- Filtration membranes

High-performance colony counters


Scan works for every kind of colony. The minimum size is **0.05 mm for Scan 1200** and **0.1 mm for Scan 300 and Scan 500**. **Scan** colony counter automatically separates confluent colonies, allows you to create polygonal exclusion areas and ignores agar flaws and air bubbles. You can also add or remove colonies manually. Every change is automatically saved in your report.



Scan read all the colonies, even the smallest




Automatic separation of confluent colonies



Automatic elimination of counting grids



Cross on each counted colony



Polygonal exclusion areas




Automated software with manual control

Color detection & chromogenic media

Scan 500 and **Scan 1200** can read chromogenic agar and differentiate colonies by color: **up to 4 different colors on the same Petri dish**. Color selection can be made directly from the color of the bacteria and a cursor allows you to set the sensitivity.

Chromogenic media reading allows the detection of *Salmonella* on XLD media and *E.Coli* on TBX media, for example.

* on Scan 500 & Scan 1200



**INCLUDED
IN THE SOFTWARE***

Scan: 3 models adapted to your needs




Scan 300 Essential

Ref 436 300

- 6 combinations of lighting and backgrounds
- Motorized background color
- Brightness, contrast and sensitivity are automatically optimized by the software
- Long lasting LED lighting
- CMOS color camera, zoom x28, M12 lens
- Minimum size of detected colony: 0.1 mm

Count these media



Surface/pour plating Spiral plating




Scan 500 Efficient

Ref 436 000

- 6 combinations of lighting and backgrounds
- Motorized background color
- Brightness, contrast and sensitivity are automatically optimized by the software
- Long lasting LED lighting
- CMOS color camera, zoom x28, M12 lens
- Minimum size of detected colony: 0.1 mm
- **Detects and counts up to 7 colors on the same dish**
- **Inhibition zone reading with EUCAST, CA-SFM, CLSI and editable database**

Count these media



Surface/pour plating Spiral plating
Chromogenic Petri dishes Antibiogram




Scan 1200 High-Resolution

Ref 437 000

- 6 combinations of lighting and backgrounds
- Motorized background color
- Brightness, contrast and sensitivity are automatically optimized by the software
- Long lasting LED lighting
- **Color HD CCD camera, zoom x28, HD Japanese lens**
- **Minimum size of detected colony: 0.05 mm**
- **Detects and counts up to 7 colors on the same dish**
- **Inhibition zone reading with EUCAST, CA-SFM, CLSI and editable database**
- **Petrifilm™, Compact Dry™, MC-Media Pads™, Easy Plate™ and filtration membrane reading**

Count these media



Surface/pour plating Spiral plating
Chromogenic Petri dishes Antibiogram
Petrifilm™ Filtration membranes Circle
Compact Dry™ MC-Media Pads™ Easy Plate™


Inhibition zone from paper discs, agar wells & peni cylinders

Performance and flexibility

Scan 500 and **Scan 1200** allow efficient work flow because you can create and edit a list of antibiotics, useful for routine analysis.

Measured by **Scan**, inhibition zones guarantee repeatability and reproducibility of analysis and diagnosis reliability.

- Rapid detection: up to 8 antibiotic sensitivities in 1 click.
- Paper discs, agar wells and peni cylinders may be manually added or deleted. Inhibition zones may be manually resized.



The result of sensitivity in contact with the antibiotic is fast and visualization of results is clear:

- > Red (resistant)
- > Yellow (intermediate)
- > Green (susceptible)

Medical analysis


Inhibition zone measurement allows you to test the efficiency of antibiotics on micro-organisms to accelerate the diagnosis in order to choose precisely an appropriate antibiotic treatment for a patient. **Scan** has a built-in antibiotic database from the **French Society of Microbiology (CA-SFM)**, the **European Committee on Antimicrobial Susceptibility Testing (EUCAST)** and the **Clinical and Laboratory Standards Institute (CLSI)** which determines the sensitivity of the bacteria to the antibiotic. This database is fully editable.



Precision of inhibition zone radius measurement from paper discs: 0.3 mm

Pharmaceutical analysis

In the pharmaceutical industry, **Scan** allows you to test the quality of an antibiotic during its manufacturing process by measuring the inhibition zones. To evaluate the action of an antibiotic, antibiotic diffusion from paper disc, agar well or peni cylinder is supported.



Precision of inhibition zone radius measurement from agar wells: 0.3 mm

21 CFR
Part 11



* on Scan 500 & Scan 1200

Comfort of use

> High definition live image


This feature enables total control of colony counting.

Optimal visualization


Enjoy comfortable viewing of the colonies with the unequalled **Dark Field technology**, high definition live image and with the automatic optimization of the image (lighting, contrast and sensitivity). You can also check key areas thanks to the digital zoom.



Dark Field: LED are disposed in a circle for optimal contrast



Scan automatically optimizes contrast, luminosity and sensitivity



Digital zoom with the mouse wheel (up to x28)





Image with Scan 1200

Easy-to-use

All **Scan** functions are in **one single window** and colonies are counted in one click.

The **Scan** easy commands (visualization, settings and results) allow quick access to both ongoing and archived work sessions.

Scan software is available in **7 languages** (English, French, Chinese, Russian, Japanese, Spanish and German) and is updated regularly. The intuitive use of **Scan** does not require any special training.



Fast communication, total traceability

Results harmonization

Using the **Scan** allows more reliable analyses and harmonizes the results within a team.

You can save as many settings as you wish and customize the settings according to the type of dishes and agar you use.

The automatic archiving of data, photos, comments and results ensures total traceability.



Print your results

You can export your results to your PC, archive it in Excel™, PDF, SCA or BIO format. You can also export pictures from the camera in JPEG, PNG and BMP format.

interscience

Add your own logo in the reports

Petri dish before counting

Petri dish after counting

Sample

Sample analysed with SCAN 1200®, version 7.0.2.0

CFU/mL : 1.12E+02
CFU : 96
Ø 79.6 mm

Sample information

Operator name :	be	Sample N° :	2	CFU/mL :	1.12E+02
Parameters :	Coli-Coliforms	Count :	96	Dilution :	1
Date Time :	12/10/2014 9:52:13	Area (%) :	86 %		

Comments

Comment :
OK

E. Coli : 23 CFU 2.68E+01
Coliform : 73 CFU 8.51E+01


Analysis result

Printed report example

Reproducibility of results

Automatic counting is a guarantee of **regularity and standardization** of analyses, which is the key to ensure accurate and reliable results. **Reproducibility** of results is guaranteed whatever the day, conditions and user.

A scientific study has proven that **Scan** colony counters can achieve an accuracy of up to 98%.



Manual counting:

Random results over time and different users by manual counting of colonies

Automatic counting:

Standardized and reproducible results by automatic counting

Study made on *Bacillus cereus*, *Escherichia coli* and *Lactobacillus casei*

Internal traceability


Thanks to the **LIMS** connection and the barcode reader, photos of counted plates are **saved and traceable**. The images are accessible and recountable at any time.



Barcode



Archiving




LIMS connection



dataLink
(see p. 12)

Secure your sessions

Sessions are secured with a **security code** (one per operator) and the impossibility to alter each saved counting. **Scan** use allows the compliance with **21 CFR part 11**: system securization, operational controls and documentation management.



> External traceability

Scan software provides numerous possibilities to easily and quickly export your results.




Work sessions saving




PDF export



JPEG, PNG & BMP
formats export



Export results to
Excel™ to ensure
traceability



Print report from Scan

Plate & Count System with dataLink

INCREASE
your lab capacities with
easySpiral & Scan!



Plate & Count System with dataLink enables automatic plating and colony counting with full traceability!

- **INCREDIBLE SAVINGS:** Save up to 75% in time, consumables and bench space
- **FAST:** Full plating cycle in 25 seconds and counting in 1 click. No manual data input.
- **RELIABLE:** 98% repeatable and reproducible results
- **TOTAL TRACEABILITY WITH dataLink:** Automatic data saving and reporting

Plate & Count System with dataLink includes:

- **easySpiral**: Automatic Spiral platers
- **Scan**: Automatic colony counters
- **dataLink**: Petri dish labeling solution


HOW DOES IT WORK?

STEP 1



Plating with easy**Spiral Pro** or easy**Spiral Dilute**. The supervision software collects data from the plater.

STEP 2



Label printing with the Datamatrix code. Stick the label on the plated Petri dish before incubation.

STEP 3



* Please check LIMS compatibility

... Incubation 24-72 h

After incubation, scan the Datamatrix code. The **Scan** colony counter is parametered automatically with the Datamatrix data. Click on "VALIDATE". Export the data.

PLATE AND COUNT YOUR PETRI DISHES

■ From 30 to 1×10^{12} cfu/mL on one single Petri dish

easy**Spiral** automatically plates a sample in **8 seconds**: from 30 to 1×10^{12} CFU/mL on a **single Petri dish** without prior sample dilution. Once the sample is plated and incubated, it is ready to be counted by **Scan** automatic colony counters. Results are immediately displayed and saved.

■ Up to 75% savings

easy**Spiral** and **Scan** guarantee the regularity and standardization of the analyses, save time, consumables and bench space of up to 75%.



Technical specifications

	Scan 300		Scan 500		Scan 1200	
Reference	436 300		436 000		437 000	
OVERVIEW						
Painted stainless steel shell	✓		✓		✓	
LIMS/SIL connection	✓		✓		✓	
USB connection	✓		✓		✓	
Available with dataLink/dataLink pro traceability system	✓		✓		✓	
Counting on pour, surface, Spiral and circle plated Petri dishes	✓		✓		✓	
Counting on chromogenic dishes	-		✓		✓	
Counting on Petrifilm™, Compact Dry™, MC-Media Pad™, Easy Plate™, filtration membranes	-		-		✓	
Automatic counting	✓		✓		✓	
Inhibition zone reading	-		✓		✓	
COUNTING						
Counting	Automatic with manual control (adding/removing colonies)					
Automatic separation of clustered colonies	✓		✓		✓	
Creation of polygonal exclusion zones	✓		✓		✓	
Counting time	Up to 1000 colonies per second					
Minimal size of colony	0.1 mm				0.05 mm	
INHIBITION ZONE READING						
Antibiotic disc detection	-		Automatic with possibility to manually add or remove			
Automatic detection of antibiogram support	-		Disks (several brands simultaneously), wells			
Display resolution	-		± 0.1 mm			
Inhibition zone measurement accuracy	-		± 0.3 mm			
Number of antibiotic paper disks	-		Up to 7 antibiotics on a Ø 90 mm Petri dish			
Reading time	-		7 inhibition zone reading between 1 to 3 s			
Interpretation system	-		CA-SFM Human health / EUCAST / CA-SFM Veterinary / CLSI (Clinical, Laboratory Standards Institute) / Customizable list			
SPECIFICATIONS						
Color camera	CMOS				HD CCD	
Lens	M12 Lens				HD Japanese lens	
Zoom	x 28					
Resolution (megapixels)	1				1.2	
White LED Lighting technology	Dark Field					
LED Lighting system	Automatic with 6 combinations, top and/or bottom light, white or black background					
Petri dish size	Ø 55 mm - Ø 90 mm					
Color detection	-		4 colors on the same dish + 2 colors to exclude			
Voltage - Frequencies	100-240 V~ 50-60 Hz					
Warranty	3-year (after recording the warranty)					
Spare parts availability	10 years					
In compliance with	21 CFR Part 11, ISO 7218 and AOAC 977.27					
SOFTWARE						
Languages	English, French, Japanese, Chinese, Russian, Spanish, German					
License type	Perpetual					
Software updates	3 years excluding qualification (after recording the warranty)					
User account management	Local or LDAP					
TRACEABILITY						
USB Data export	Recountable session, Excel™, PDF report, jpg, png and bmp images					
Data security	Modified data traceability in conformity with 21 CFR part 11					
Results/traceability	Image / sample number / comments /date / time		Image / sample number / comments / date / time / antibiotic name / bacterial name / measured diameter / result categorized according to standards / minimum and maximum critical diameter			
PC MINIMUM REQUIREMENTS						
Operating sytem	Windows™ 10 or 11 (or higher)					
Processor	Intel i5, 13th or 14th generation minimum					
RAM	4 GB for use of the Scan					
Equipments	USB port free					
Screen	1280 X 1024 pixels or higher					

Scan 300: delivered with 1 Scan software, 1 15V power supply, 1 power supply, 1 USB cable, 2 validation plates, 1 user manual, 1 conformity certificate / warranty card.

Scan 500 / Scan 1200: 1 Scan software, 1 15V power supply, 1 power supply, 1 USB cable, 3 validation plates, 1 user manual, 1 conformity certificate / warranty card.

PC requirements are subject to change. Please check our website www.interscience.com for current updates and additional informations.

Certified production



Product made for INTERSCIENCE by Interlab, an ISO 9001 certified company.

Scan accessories



Barcode reader
Ref.: 522 000




**Adaptor for
Petri dish (55 mm)**
Ref.: 436 005




**Adaptor for
MC-Media Pads^{TM***}**
Ref.: 437 001



**Adaptor for
Petrifilm^{TM***}**
Ref.: 437 002



**Adaptor for
Easy Plate^{TM***}**
Ref.: 437 003



**Adaptor for
Compact Dry^{TM***}**
Ref.: 437 004



**Validation plates
for counter and
reader^{**}**
Ref.: 437 006




**Certified validation
plates for counter
and reader^{**}**
Ref.: 437 007



**Validation plates
for Scan 300^{*}**
Ref.: 437 008

* On Scan 300 / ** On Scan 500 & Scan 1200 / *** On Scan 1200

Discover our complete range for microbiology



www.interscience.com



www.interscience.com



interscience

PARIS

Phone: +33 (0)1 34 62 62 61
Email: info@interscience.com

FRANKFURT

Phone: +49 611 7238 7770
Email: sales.germany@interscience.com

BOSTON

Phone: +1 781 937 0007
Email: sales.usa@intersciencelab.com

SHANGHAI

Phone: +86 (0)21-64739390
Email: sales.china@interscience.cn

SINGAPORE

Phone: +65 6977 7232
Email: sales.asia@interscience.com

TOKYO

Phone: +81 3 6712 9715
Email: sales.japan@interscience.com

Your local distributor