



3 dilutions on only 1 Petri dish!

interscience













. . . .

1979

Spiral® DS

1996

Spiral® DS+

2010

easy**Spiral**® easy**Spiral Pro**®

2012

easySpiral Dilute®

2020

easySpiral Pro® Milk


45 years of **Spiral®** method experience


- 1973 The Spiral® method is invented by Dr. Ed Campbell, researcher at the FDA (Food and Drug Administration) in Cincinnati, USA.
- **1992** Dr Ed. Campbell and his friend François Jalenques, founder of intenscience, patent together an updated version of the Spiral® method in France.
- **2010** interscience introduces the easy**Spiral**® and easy**Spiral** Pro® with a brand-new architecture: they are faster, simpler and cheaper than the competition.
- 2012 Launching of easySpiral Dilute® is the 1st automatic diluter and plater: quickly becomes our best-selling Spiral® plater interscience.
- **2019** The introduction of the easy**Spiral Pro® Milk** expands the widest range on the market.

Since 45 years, intenscience's automatic Spiral® platers have been a reference for food, medical, public institutes research, cosmetics, pharma and probiotics labs throughout the world. ISO 7218, ISO 4833-2, FDA-BAM and AOAC 977.27 standards recommends the use of Spiral® platers for bacterial analysis.

What is Spiral® method?

With easySpiral® automatic platers, increase your lab capacity with a counting range from 100 to 1x10⁷ CFU/mL on 1 single Petri dish without prior dilution.





How to count?

The sample is deposited on the surface of a rotating Petri dish by a precise measuring microsyringe, according to a **logarithmically decreasing Archimedes spiral** from the center of the dish to its periphery. The volume is calibrated and known at every point of the Petri dish. Bacterial concentration is determined by dividing the number of colonies found by the volume dispensed in the same sector of the dish, either manually, or by using an automatic colony counter.

Manual colony counting

Scan® 100 with Spiral® counting grid

Automatic colony counting

Scan® 300 | 500 | 1200 | 4000

Real-time incubation and colony counting

ScanStation® 100 | 200 | 300







Key features

Compact

40 cm wide, can be used inside and outside a laminar flow Stainless steel housing

Automatic disinfection

Overflow technology cleaning system No cross contamination

PATENTED

Easy serial dilutions*

Choose your dilution factor

Do 10⁻⁵ dilution or more automatically

No more 9 mL tube preparation

PATENTED

One-touch plating

Intuitive control panel
Choose from 4 plating modes*
Plate on Ø 55, 90 and 150 mm*
Petri dishes

USB programmable volume & traceability*

Easy PC connectivity via USB
From 10 to 1000 microliters
programmable volumes
3 preset plating modes
Monitoring software with traceability


COMPLIANT WITH













Direct sampling from a tube on easySpiral Pro® Milk*

Ideal for raw materials and the milk payment Direct intake from the holder/pill container



Glass syringe typical accuracy of 0.5% No more pipet calibration



2 fully autoclavable bottles and CPC connectors Alarm when bottle is empty

High speed automation

Fast rotating arm
Plating cycle: 25 s
1 disinfection, 5 dilutions
and 1 plating in 134 s*

PATENTED










	easy Spiral ®	easy Spiral Pro ®	easy Spiral Pro® Milk	easy Spiral Dilute®
Full plating cycle in 25 s	✓	✓	✓	✓
Controlled by microprocessor	✓	✓	✓	✓
Overflow stylus disinfection system (patented)	✓	✓	✓	✓
Programmable dispensed volume via USB	-	✓	✓	✓
Traceability	-	✓	✓	✓
Compatible with dataLink™	-	✓	✓	✓
Sampling from a tube	-	-	✓	-
Dilution mode (patented)	-	-	-	✓
Counting range up to 1x10 ¹² CFU/mL	-	-	-	✓
Sample mixing before dilution/plating	-	-	-	✓

easy**Spiral** Dilute®

World 1st and only automatic diluter & plater



easy**Spiral** Dilute® is a **2-in-1 automatic diluter and plater**, which allows you to do **5 x 1/10th serial dilutions**. It will then automatically plate on 1 single Petri dish, with a countable range from **30 to 1 x 10¹² countable CFU/mL**. **No more manual dilutions and plating is automatic!**

- ▶ Place your sample in the beaker
- 1 The easy Spiral Dilute® will automatically dilute the sample up to 5 times
- 2 The easy Spiral Dilute® will automatically plate the sample on a Petri dish from 30 to 1x1012 countable CFU/mL

Questions

Which diluents can I use?

easy**Spiral** Dilute® works with saline solutions (e.g.: Ringer's solution, physiological saline).

Avoid connecting diluents with nutrients to the easy**Spiral** Dilute® as they may cause contaminations. Bacteria only remain a few seconds in the beakers before being deposited on the Petri dish so using a saline diluent has no effect on the final enumeration.

It is however possible to use nutrient solutions as diluent in the semiautomatic "My Diluent" mode. In that case, fill in the beakers with 4.5 mL diluent before processing.

I need to make a 10⁻⁷ dilution. How can I do it?

The easy Spiral Dilute® can do as many dilutions as you wish.

Proceed with a 10^{-5} dilution, take the last beaker and place it as the mother sample and then proceed with a new automatic dilution sequence...

The accuracy of the volumes plated, intaken or distributed by the easySpiral Dilute® is high: 2% error. How can I check it?


It is easy and automatic.

Test of the deposited volume during plating:

Press on the plating volume test button. This function will fill a previously weighed beaker with 500 μ L. Weigh the distributed volume to check the results.

Test of the volume of the dilutions:

Press on the dilution volumes test button. This function will fill 2 previously weighed beakers: one with 0.5 mL and the 2nd one with 4.5 mL. Weigh the distributed volumes to check the results.




Why use an automatic diluter & plater?

easy**Spiral** Dilute® allows you to save 50 % on consumables and 50 % on your time.

Comparison between manual and automatic method below. Example of 50 samples: Dilution: 10^{-7}

Plating: 10⁻⁵, 10⁻⁶ and 10⁻⁷



Plating modes

Exponential mode

Decreasing surface concentration

 $50~\mu L,\,100~\mu L,\,200~\mu L$ and USB programmable

Counting on Petri dish from 30 to 1x10¹² CFU/mL

Constant⁽¹⁾

Homogeneous surface deposition

 $50~\mu\text{L},\,100~\mu\text{L},\,200~\mu\text{L}$ and USB programmable

Easy colony counting: No colonies on the edge of the Petri dish

Circle mode

3 or 6 dilutions on 1 Petri dish

10 μL inner circle, 16 μL middle circle and 25 μL outer circle


Savings on Petri dishes

Pour⁽¹⁾ mode

Homogeneous deposition in a plate without agar⁽³⁾

 $50~\mu L,\,100~\mu L,\,200~\mu L$ and USB programmable

Easy colony counting for invasive colonies. Pharmacopoeia compliant



⁽¹⁾ Depending on model (see technical specifications p11) (3) In pour mode, easySpiral Dilute® deposits one sample drop at the bottom of an empty Petri dish. Then, the user adds the molten agar on top and homogenizes it with the sample.

Choose your automatic plater








			easy Spiral °	easy Spiral Pro ° & easy Spiral Pro ° Milk	easy Spiral Dilute®
			Counting range from: 300 to 1.3x10 ⁵ CFU/mL	Counting range from: 100 to 1x10 ⁷ CFU/mL	Counting range from: 30 to 1x10 ¹² CFU/mL
.es		w stylus disinfection (patented)	✓	✓	✓
Top features	Traceak	bility	-	✓	✓
	Dilution mode (patented)		-	-	✓
	0	Exponential plating mode	✓	✓	✓
Plating modes		Circle plating mode (patented)	✓	✓	✓
		Constant plating mode	-	✓	✓
	•	Pour plating mode	-	-	✓
size	0	90 mm Petri dishes	✓	✓	✓
Petri dish size	0	150 mm Petri dishes	-	✓	✓
	0	55 mm Petri dishes	-	-	✓

easy**Spiral**® technology

Delivering the fastest cycle time ever, easySpiral® has an impressive rotating arm architecture, patented disinfection system avoiding risks of cross-contamination and full traceability of the operations with a counting range from 30 to 1x10¹² CFU/mL on 1 Petri dish.





Only 25 seconds for a plating cycle!

- 25 seconds for a plating cycle (disinfection, sample-taking and plating)
- Rotating arm with high speed movement
- Patented electronic robot system
- 10 Petri dishes plated in less than 2 mn (with the same sample)

Overflow technology®

- No cross-contamination
- Unique patented disinfection system by overflow
- 600 cleanings without changing bottles (2L)





Accessible fluidics

- Change the stylus in 30 seconds
- Visible and accessible connections in the event of a leak
- Removable bathtub

Traceability with dataLinkTM

 All the sample data are set on the Datamatrix label (volume, dilution, plating mode, lighting settings...).



Technical specifications

		easy Spiral ®	easy Spiral Pro ®	easy Spiral Pro® Milk	easy Spiral Dilute®		
	Reference	412 000	413 000	413 090	414 000		
OVERVIEW	Number of dilutions at 1/10th for 1 cycle	-	-	-	5		
	Syringe capacity	1 000 μL					
	Autonomy in disinfection (in normal mode)	1000 (with 2 L bottles)					
	Autonomy in distilled water (in normal mode)	600 cycles (2L bottles) -					
	Autonomy in diluents	-	-	-	63 cycles (2 L bottle)		
	Approved diluents	-	-	-	Saline water Ringer solution		
	Pressure in the stylus	up to 8 bars					
	Traceability	- Excel™, LIMS, dataLink™					
PLATINGS	Exponential mode (Spiral® mode): Petri dishes diameters (Ø) available	90 mm	90 mm a	nd 150 mm	90 mm and 150 mm		
	Circle mode: Petri dishes diameters (Ø) available	90 mm	90 mm and 150 mm		3 dilutions on 90 mm and 6 dilutions on 150 mm		
	Constant mode: Petri dishes diameters (Ø) available	-	90 mm and 150 mm		55 mm, 90 mm and 150 mm		
	Pour plate mode: Petri dishes diameters (Ø) available	-	-	-	55 mm, 90 mm and 150 mm		
	Counting range	300 to 1.3x105 UFC/mL	100 to 1x	100 to 1x10 ⁷ UFC/mL 30 to 1x10 ¹² U			
	Programmable dispensed volume via USB	- from 10 μL to 1 000 μL					
	Preset dispensed volume	50 μL and 100 μL	100 μL and 200 μL				
	Successive plating capacity with the same sample	up to 20 Petri dishes (50 μL)					
	Sample intake + 1 disinfection + 1 plating	25 s					
TIME	1 disinfection + 5 dilutions + 1 plating of 1 dilution	-	-	-	134 s		
	1 disinfection + 5 dilutions + 1 plating of each dilution	-	-	-	234 s		
	Controlled by microprocessor	✓	✓	✓	✓		
	All stainless steel	✓	✓	✓	✓		
S	1 year warranty	✓	✓	✓	✓		
읽	Disinfection process	✓	✓	✓	✓		
5	Filling time and volume programmable via USB	-	✓	✓	✓		
SPECIFICATIONS	Sampling from beakers (DB50)	✓	✓	-	✓		
ß	Sampling from tubes (Dim: 80 x 35 mm)	-	-	✓	-		
	My Diluent key	-	-	-	✓		
	Sample mixing before dilution/plating	-	-	-	✓		
28	Dimensions (w x d x h)	40 x 41.5 x 29 cm					
SIO	Weight	16.4 kg					
WEIGHT AND DIMENSIONS	Voltage - frequencies	100-240V~ 50-60 Hz					
₹₫	Power	45 W					

easySpiral Pro*, easySpiral Pro* easySpiral Pro* Milk, easySpiral Dilute* are delivered with 1 turntable ring for Ø 90 mm Petri dishes (by default on the unit), 1 blue dye for testing, 1 syrings with silicone tubing, 1 plating ring for Ø 55 mm Petri dishes*, 1 plating ring for Ø 150 mm Petri dishes*, 1 Spiral* counting grid for Ø 90 mm Petri dishes*, 1 circle mode counting grid for Ø 90 mm Petri dishes*, 1 circle mode counting grid for Ø 90 mm Petri dishes*, 1 circle mode counting grid for Ø 150 mm Petri dishes*, 1 circle mode c

Certified production









Product made for INTERSCIENCE by Interlab, an ISO 9001 certified company.


Accessories







for Spiral® platers (5 mL) Ref 415 100



⁽¹⁾ Depending on model (see technical file) (2) Not available on easy Spiral Pro® Milk



The Plate & Count system® with automatic diluting, plating and colony counting are now the perfect solution when it comes to meeting efficiency and traceability requirements for microbiological analyses.

- 75% time-saving and consumables savings
- Secured and traceable bacterial results
- Top-of-the-line technology made in France

For more info on Plate & Count system, please refer to the brochure.

Your local distributor

interscience in the world

Phone: +33 (0)1 34 62 62 61 - Email: info@interscience.com

FRANKFURT

Phone: +49 611 7238 7770 - Email: sales.germany@interscience.com

Phone: +1 781 937 0007 - Email: sales.usa@intersciencelab.com

SHANGHAI

电话: +86 (0)21-64739390 - 邮址: sales.china@interscience.com

SINGAPORE

Phone: +65 6977 7232 - Email: sales.asia@interscience.com

