

Guide des tubulures

A Ensemble de distribution

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E : Excellent / G : Good / F : Fair / U : Not Recommended

Environment, % Conc.* w-Water alc-Alcohol	Tygon™ 3350	Silicone	Environment, % Conc.* w-Water alc-Alcohol	Tygon™ 3350	Silicone	Environment, % Conc.* w-Water alc-Alcohol	Tygon™ 3350	Silicone
Acetaldehyde	F	F	Calcium Salts	E	E	Gallic Acid, 17% in acetone	U	U
Acetamide, 67% in w	F	G	Calcium Sulfate, 1% in w	E	E	Gasoline, Automotive	U	U
Acetate Solvents	U	U	Carbon Dioxide, Wet/Dry	E	E	Gelatin	E	E
Acetic Acid, 10% in w	E	E	Carbon Disulfide	U	U	Glucose, 50% in w	E	E
Acetic Acid, 50-60% in w	E	E	Carbon Monoxide	E	E	Glycerol, (Glycerin)	E	E
Acetic Acid, Glacial, 100%	U	U	Carbon Tetrachloride	U	U	Glycolic Acid, 70% in w	E	E
Acetic Anhydride	E	F	Carbonic Acid	E	G	Heptane	U	U
Acetone	F	U	Castor Oil	E	E	Hexane	U	U
Acetonitrile	U	U	Cellosolve	U	U	Hydrazine	U	U
Acetyl Bromide	U	U	Cellosolve Acetate	U	U	Hydrobromic Acid, 20-50% in w	U	U
Acetyl Chloride	U	U	Chlorine, Dry Gas	U	U	Hydrobromic Acid, 100% in w	U	U
Acetylene Gas	E	F	Chlorine, Wet Gas	U	U	Hydrochloric Acid, 10% in w	E	E
Acrylonitrile	U	U	Chloroacetic Acid, 20% in w	E	G	Hydrochloric Acid, 37% in w	U	U
Adipic Acid, 100% in alc	U	U	Chlorobenzene, Mono, Di, Tri	U	U	Hydrocyanic Acid	E	E
Air	E	E	Chloroform	U	U	Hydrofluoric Acid, 10% in w	U	U
Alcohols General	G	E	Chlorosulfonic Acid	U	U	Hydrofluoric Acid, 25% in w	U	U
Aliphatic Hydrocarbons	U	U	Chromic Acid, 10-20% in w	U	U	Hydrofluoric Acid, 40-48% in w	U	U
Allyl Alcohol	U	U	Chromic Acid, 50% in w	U	U	Hydrogen Gas	E	E
Alum, 5% in w	E	E	Citric Acid, 10-20% in w	E	E	Hydrogen Peroxide, 3% in w	E	E
Aluminum Chloride, 53% in w	E	E	Coconut Oil	E	E	Hydrogen Peroxide, 10% in w	E	E
Aluminum Hydroxide, 2% in w	E	E	Corn Syrup	E	E	Hydrogen Peroxide, 30% in w	E	E
Aluminum Sulfate, 50% in w	E	E	Cottonseed Oil	E	E	Hydrogen Peroxide, 90% in w	F	F
Aluminum Salts	E	E	Cresol (m, o, or p)	G	G	Hydrogen Sulfide	E	E
Amines	U	U	Cresylic Acid	U	U	Hydroquinone, 7% in w	G	F
Ammonia Gas	U	U	Cupric Chloride, 40% in w	E	E	Hypochlorous Acid, 25% in w	E	E
Ammonia, Anhydrous Liquid	U	U	Cupric Nitrate, 70% in w	E	E	Iodine, 50 ppm in w	E	E
Ammonium Acetate, 45% in w	E	E	Cupric Sulfate, 13% in w	E	E	Isobutyl Alcohol	U	U
Ammonium Carbonate, 50% in w	E	E	Cyclohexane	U	U	Isocetane	U	U
Ammonium Hydroxide, 5-10% in w	U	U	Cyclohexanone	U	U	Isopropyl Acetate	U	U
Ammonium Hydroxide, 30% in w	U	U	Detergent Solutions	E	G	Isopropyl Alcohol	U	U
Ammonium Persulfate, 30% in w	E	E	Dibutyl Phthalate	E	E	Isopropyl Ether	U	U
Ammonium Salts	E	E	Diesel Fuel	U	U	Jet Fuel, JP8	U	U
Ammonium Sulfate, 30% in w	E	E	Diethylamine, 2.5% in w	U	U	Kerosene	U	U
Amyl Acetate	U	U	Diethylene Glycol	E	E	Ketones	U	U
Amyl Alcohol	U	U	Dimethylformamide	E	G	Lacquer Solvents	U	U
Amyl Chloride	U	U	Dimethylsulfoxide	F	U	Lactic Acid, 3-10% in w	E	E
Aniline	U	U	Diethyl Phthalate	E	E	Lactic Acid, 85% in w	U	U
Aniline Hydrochloride	U	U	Dioxane	U	U	Lard, Animal Fat	E	E
Antimony Salts	E	E	Ether	U	U	Lead Acetate, 35% in w	E	E
Aqua Regia	U	U	Ethyl Acetate	U	U	Lead Salts	E	E
Aromatic Hydrocarbons	U	U	Ethyl Alcohol (Ethanol)	G	F	Lemon Oil	U	U
Arsenic Acid, 20% in w	F	U	Ethyl Benzoate	U	U	Limonene-D	U	U
Arsenic Salts	E	E	Ethyl Chloride	U	U	Linoleic Acid	G	F
ASTM Reference No. 1 Oil	E	G	Ethyl Ether	U	U	Linseed Oil	E	E
ASTM Reference No. 2 Oil	G	G	Ethylene Bromide	E	E	Lubricating Oils, Petroleum	G	G
ASTM Reference No. 3 Oil	U	U	Ethylene Chlorohydrin	G	G	Magnesium Carbonate, 1% in w	E	E
Barium Carbonate, 1% in w	E	E	Ethylene Dichloride	U	U	Magnesium Chloride, 35% in w	E	E
Barium Hydroxide, 5% in w	E	E	Ethylene Glycol	E	E	MagnesiumHydroxide,10%indil.acids	E	E
Beer	E	E	Ethylene Oxide	E	F	Magnesium Nitrate, 50% in w	E	E
Benzaldehyde	F	F	Fatty Acids	G	F	Magnesium Sulfate, 25% in w	E	E
Benzene	U	U	Ferric Chloride, 43% in w	E	E	Maleic Acid, 30% in w	G	F
Benzenesulfonic Acid	U	U	Ferric Nitrate, 60% in w	E	E	Malic Acid, 36% in w	E	E
Benzoic Acid	U	U	Ferric Sulfate, 5% in w	E	E	Manganese Salts	E	E
Benzyl Alcohol	E	E	Ferrous Chloride, 40% in w	E	E	Mercuric Chloride, 6% in w	E	E
Bleach Liquor, 22% in w	U	G	Ferrous Sulfate, 5% in w	E	E	Mercuric Cyanide, 8% in w	E	E
Borax, 6% in w	E	G	Fluoboric Acid, 48% in w	U	U	Mercury	E	E
Boric Acid, 4% in w	E	G	Fluorine Gas	U	U	Mercury Salts	E	E
Bromine, Anhydrous Liquid	U	U	Fluosilicic Acid, 25% in w	F	U	Methane Gas	E	F
Butadiene	E	E	Formaldehyde, 37% in w	F	F	Methyl Alcohol (Methanol)	G	G
Butane	E	E	Formic Acid, 25% in w	E	G	Methyl Acetate	U	U
Butyl Acetate	U	U	Formic Acid, 40-50% in w	E	F	Methyl Bromide	U	U
Butyl Alcohol	U	U	Formic Acid, 98% in w	E	F	Methyl Chloride	U	U
Butyric Acid	U	U	Freon 11	E	E	Methyl Ethyl Ketone	U	U
CalciumCarbonate,25%indiluteacids	E	E	Freon 12	E	E	Methyl Isobutyl Ketone	U	U
Calcium Chloride, 30% in w	E	E	Freon 22	E	E	Methylene Chloride	U	U
CalciumHydroxide,10%inglycerol	E	G	Fruit Juice	E	E	Methyl Methacrylate	U	U
Calcium Hypochlorite, 20% in w	U	G	Fuel Oil	U	U	Milk	E	E
Calcium Nitrate, 55% in w	E	E	Furfural	U	U	Mineral Oil	U	U

Environment, % Conc.* w-Water alc-Alcohol	Tygon™ 3350	Silicone
Mineral Spirits	U	U
Molasses	E	E
Monoethanolamine	U	U
Motor Oil	U	U
Naphtha	U	U
Naphthalene	U	U
Natural Gas	E	F
Nickel Chloride, 40% in w	E	E
Nickel Nitrate, 75% in w	E	E
Nickel Salts	E	E
Nickel Sulfate, 25% in w	E	E
Nitric Acid, 10% in w	F	F
Nitric Acid, 35% in w	U	U
Nitric Acid, 68-71% in w	U	U
Nitrobenzene	U	U
Nitromethane	U	U
Nitrous Acid, 10% in w	G	G
Nitrous Oxide	E	G
Oils, Animal	E	E
Oils, Essential	U	U
Oils, Hydraulic (Phosphate Ester)	U	U
Oils, Hydrocarbon	G	G
Oils, Vegetable	E	E
Oleic Acid	G	F
Oleum, 25% in w	G	G
Ortho Dichlorobenzene	U	U
Oxalic Acid, 12% in w	E	F
Oxygen	E	E
Ozone, 300pphm	E	U
Palmitic Acid, 100% in ether	G	F
Paraffins	U	U
Perchloric Acid, 67% in w	U	U
Perchloroethylene	U	U
Phenol, 5-10% in w	E	U
Phenol, 91% in w	G	U
Phosphoric Acid, <10% in w	F	U
Phosphoric Acid, 25% in w	U	U
Phosphoric Acid, 85% in w	U	U

Environment, % Conc.* w-Water alc-Alcohol	Tygon™ 3350	Silicone
Phosphorous Trichloride Acid	U	U
Photographic Solutions	G	F
Phthalic Acid, 9% in alc	G	F
Phthalic Anhydride, 9% in alc	E	F
Picric Acid, 1% in w	U	U
Plating Solutions	U	U
Potassium Carbonate, 55% in w	E	E
Potassium Cyanide, 33% in w	E	E
Potassium Dichromate, 5% in w	E	E
Potassium Hydroxide, <10% in w	E	G
Potassium Iodide, 56% in w	E	E
Potassium Permanganate, 6% in w	E	E
Potassium Salts	E	E
Propane Gas	E	F
Propylene Glycol	E	E
Propylene Oxide	E	E
Pyridine	U	U
Salicylic Acid, 1% in w	E	E
Silicone Oils	U	U
Silver Nitrate, 55% in w	E	E
Skydrol 500A	U	U
Soap Solutions	E	F
Sodium Acetate, 55% in w	U	U
Sodium Benzoate, 22% in w	E	E
Sodium Bicarbonate, 7% in w	E	E
Sodium Carbonate, 7% in w	E	E
Sodium Chlorate, 45% in w	E	E
Sodium Chloride, 20% in w	E	E
Sodium Cyanide, 30% in w	E	G
Sodium Fluoride, 3% in w	E	E
Sodium Hydroxide, 10-15% in w	E	E
Sodium Hydroxide, 30-40% in w	E	G
Sodium Hypochlorite, 5.5% in w	U	E
Sodium Hypochlorite, 12.2% in w	U	G
Sodium Nitrate, 3.5% in w	E	E
Sodium Salts	E	E
Sodium Sulfate, 5% in w	E	E
Sodium Sulfide, 45% in w	E	E

Environment, % Conc.* w-Water alc-Alcohol	Tygon™ 3350	Silicone
Stannic Chloride, 50% in w	E	E
Stannous Chloride, 45% in w	E	E
Stearic Acid, 5% in alc	G	F
Styrene Monomer	U	U
Sulfur Chloride	U	U
Sulfur Dioxide, Gas Dry	E	G
Sulfur Dioxide, Gas Wet	E	G
Sulfur Trioxide, Wet	G	F
Sulfuric Acid, 10% in w	E	E
Sulfuric Acid, 30% in w	G	G
Sulfuric Acid, 95-98% in w	U	U
Sulfurous Acid	E	E
Tannic Acid, 75% in w	E	F
Tartaric Acid, 56% in w	E	G
Tetrahydrofuran	U	U
Thionyl Chloride	E	G
Tin Salts	E	E
Titanium Salts	E	E
Toluene	U	U
Trichloroacetic Acid, 90% in w	E	G
Trichloroethane	U	U
Triethanolamine	U	U
Trichloroethylene	U	U
Trichloropropane	U	U
Tricresyl Phosphate	E	E
Trisodium Phosphate	E	E
Turpentine	U	U
Urea, 20% in w	E	G
Uric Acid	E	G
Vinegar	E	E
Vinyl Acetate	U	U
Water, Deionized	E	E
Water, Distilled	E	E
Xylene	U	U
Zinc Chloride, 80% in w	E	E
Zinc Salts	E	E

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