


BagLight® PolySilk®

Blender bag without filter, transparent

Ref. 131 025 (100 mL) – **Ref. 132 025** (400 mL) – **Ref. 132 050** (400 mL) – **Ref. 132 200** (2000 mL)
Ref. 133 025 (3500 mL) - **Ref. 132 525 unmarked** (400 mL)



Standards:



Manufactured under:



BagLight® PolySilk® is a blender bag without filter. It is ideal for sample preparation.

TECHNICAL SPECIFICATIONS

- Non-filter bag
- PolySilk®, polyolefin complex
- Flexible and transparent
- Resistant to freezing and high temperatures (from -40 °C to 80 °C / -40 °F to 176 °F)
- No contact between the sample and the blender during homogenization of the sample
- Compatible with any lab blender

- Gamma ray treated: Gamma 5 to 12 kGy, with certificate
- Shelf-life: 10 years
- Storage: in a temperate and dry place, away from humidity and light
- Available in 100/400/2000/3500 mL

- In compliance with ISO 7218, ISO 6887 and FDA BAM (Bacteriological Analytical Manual)
- Approved for food contact - Regulation (EC) N° 1935/2004
- Designed and made in France

BagLight® PolySilk® 100

Ref. 131 025

- Max blending volume: 100 mL
- Optimal blending volume: 5 - 50 mL
- Bag dimensions: 110 x 200 mm
- Box dimensions: 31 x 22 x 7 cm, weight: 1.80 kg
- Pack of 25
- Box of 500

BagLight® PolySilk® 400

Ref. 132 025

- Max blending volume: 400 mL
- Optimal blending volume: 50 - 300 mL
- Bag dimensions: 175 x 300 mm
- Box dimensions: 35 x 21 x 11 cm, weight: 4.10 kg
- Pack of 25
- Box of 500

BagLight® PolySilk® 400

Ref. 132 050

- Max blending volume: 400 mL
- Optimal blending volume: 50 - 300 mL
- Bag dimensions: 175 x 300 mm
- Box dimensions: 35 x 21 x 11 cm, weight: 4.10 kg
- Pack of 50
- Box of 500

BagLight® PolySilk® 2000

Ref. 132 200

- Max blending volume: 2000 mL
- Optimal blending volume: 400 - 1500 mL
- Bag dimensions: 250 x 380 mm
- Box dimensions: 48 x 29 x 14 cm, weight: 10.2 kg
- Pack of 25
- Box of 500

BagLight® PolySilk® 3500

Ref. 133 025

- Max blending volume: 3750 mL
- Optimal blending volume: 400 - 3750 mL
- Bag dimensions: 380 x 560 mm
- Box dimensions: 41 x 34 x 12 cm, weight: 7.70 kg
- Pack of 25
- Box of 250

BagLight® PolySilk® 400 unmarked

Ref. 132 525

- Max blending volume: 400 mL
- Optimal blending volume: 50 - 300 mL
- Bag dimensions: 175 x 300 mm
- Box dimensions: 35 x 21 x 11 cm, weight: 4.10 kg
- Pack of 25
- Box of 500