

BagTools®

Essential accessories

A complete range of accessories for a safe and efficient laboratory work.




BagRack®

Bag holder

BagRack® allows a clean storage of your bags. You can let the bags open for easy pipeting or close it with BagClip®.

3 models: 100 ml, 400 ml and 3500 ml



Option: DripTray® 400 security drip tray in case of liquid leak, for hygiene and cleanliness.



BagOpen®

Bag opener/holder

BagOpen® allows an easy opening of the bag without contact with the sample. Ideal tool for laboratory manipulations.

3 models: 100 ml, 400 ml et 3500 ml




Thanks to BagOpen®, pipeting sample is quick and easy.



BagPipet® & BagTips®

Pipet and sterile tips

BagPipet® pipets fixed chosen volumes (0.1 ou 0.9 ou 1 ml). Used with a BagTips®, one-use sterile straws, easy to eject, BagPipet® allows the pipeting of all products, even viscous, with no blocking risks.



BagPipet® is particularly adapted to viscous and food samples pipeting, after blending in a BagFilter® filter bag with a BagMixer® blender.



BagClip®

Closing clip

BagClip® close by pinching any kind of bags, allowing a perfectly tight closing thanks to its three points closing system. BagClip® is trademarked.

3 models: 100 ml, 400 ml et 3500 ml



Thanks to its quick and perfectly tight closing, with BagClip® you can store your samples in a BagRack® in a secured way.



PetriPile®

Petri dishes storage stack

PetriPile® allow you to manipulate and transport in a secured manner all your Petri dishes. PetriPile is designed for 36 dishes (16 dishes for the 150 mm PetriPile®).

4 models: Ø 55 / 65 / 90 / 150 mm

Stainless steel designed, stackable and autoclavable, PetriPile® is ideal for incubator loading, according to AFNOR NF V 08-002 standard.

BagTools® Brochure 00311-1 - BagTools®, BagRack®, BagOpen®, BagClip®, BagPipet®, BagPipet®, BagTips®, BagFilter®, BagMixer® et PetriPile® are registered trademarks. Information in this brochure are not contractually binding.




BackRack®
Storage rack
for bags



- Price
- Loan
- Documentation

**BackPipet®
& BackTips®**
Pipet and
sterile tips



- Price
- Loan
- Documentation

BagOpen®
Bag opener/
holder



- Price
- Sample
- Documentation

BagClip®
Closing clip



- Price
- Sample
- Documentation

PetriPile®
Petri dish
storage stack



- Price
- Loan
- Documentation

Please return this form to your local distributor or to interscience.

I would like to receive more info on:

- MiniMix®** blender
- JumboMix®** blender
- BagFilter®** sterile filter-bags
- BagTools®** accessories
- Gravimat®** gravimetric dilutors
- easySpiral®** automatic platers
- Scan®** colony counters

Name

Company

Address

.....

ZIP/State/CityCountry.....

Phone Fax

Email

Web.....



www.interscience.fr

International : 30. Ch Bois Arpents, 78860 St NOM FRANCE
Tel. +33(0)1 34 62 62 61 / Fax +33(0)1 34 62 43 03
info@interscience.fr / www.interscience.fr

USA/Canada: 199 Weymouth St. ROCKLAND, MA 02370 USA
Tel. +1 781 792 2133 / Fax +1 781 792 2134
info@intersciencelab.com / www.intersciencelab.com