

## JumboMix® 3500 P CC

3500 mL lab blender

Ref. 031 110



Standards:



Manufactured under:



**JumboMix® 3500 P CC** allows blending of large samples, up to 375 g. Adjustable blending power (patented) improves bacterial extraction efficiency by adapting to the matrix. **JumboMix® 3500 P CC** is equipped with the patented Click & Clean® click-on paddle system, which facilitates cleaning by providing full access to the mixing chamber.

### TECHNICAL SPECIFICATIONS

- Stainless steel (304L) outer shell
- Adjustable blending power (Patented)
- Multi-function digital screen
- Adjustable paddles
- Side-by-side paddle stop
- Click & Clean® removable paddles (patented)
- Plain door
- 270° opening door
  
- Optimal blending volume: 200 - 3750 mL
- Q-Tight® closing: Special closing to prevent leakage from the top of the bag
- Brushless motor
- Variable speed (1.5 - 3 - 4.5 - 6 strokes/s)
- Variable blending time (1 s - 1 h or ∞ / count-down)
  
- Voltage and frequency: 100-240 V~ 50-60 Hz / 8A
- Fuses : T8A L 250V 5x20mm
- Max power: 1000 W
- International protection marking: IP21
  
- 3-year warranty (after recording the warranty card)
- Shock absorbers lifetime warranty
- Availability of spare parts : 10 years
  
- In compliance with: ISO 7218, ISO 6887-1 and FDA BAM (Bacteriological Analytical Manual)
- Manufactured under CE, RoHS, WEEE
- GLP-compliant use (Good Laboratory Practice)
- Designed and made in France

## WEIGHT AND DIMENSIONS

- Dimensions (w x d x h): 51.2 x 52.3 x 46.8 cm, net weight: 52.2 kg
- Box (w x d x h): 61 x 57 x 57 cm, weight: 56.9 kg

## DELIVERED WITH

- 1 power cable
- 1 user's manual
- 1 free pack of blender bags
- 1 certificate of conformity / warranty card

## ACCESSORIES

- Ref 036 000 **DripTray:** Stainless steel security drip tray for **JumboMix®**
- Ref 037 000 **JumboStand:** Stainless steel ergonomic table for **JumboMix® 3500**.  
Dim. (w x d x h): 58 x 51 x 54 cm. Loading height of the **JumboMix®**: 45 cm