

JumboMix® 3500 W CC

3500 mL lab blender

Ref. 032 110

WINDOW DOOR
LIFETIME
WARRANTY

SHOCK ABSORBERS
LIFETIME
WARRANTY

3
YEAR
WARRANTY



Standards:

ISO
7218

ISO
6887-1

FDA BAM
Bacteriological Analytical
Manual

Manufactured under:

RoHS
2002/95/EC

WEEE
2002/96/EC

CE UK
CA



JumboMix® 3500 W CC allows to blend large samples up to 375 g. The adjustable blending power (patented) improves the efficiency of bacterial extraction by adapting to the matrix. **JumboMix® 3500 W CC** is fitted with removable Click & Clean® paddle system (patented) that simplifies cleaning by providing full access to the blending room. The glass door allows quick controls during the blending process.

TECHNICAL SPECIFICATIONS

- Stainless steel (304 L) outer shell
- Adjustable blending power (patented)
- Multi-function digital screen
- Adjustable paddles
- Side-by-side paddle stop
- Click & Clean click-on paddles (Patented)
- Window door
- 270° opening door
- Optimal blending volume: 200 - 3750 mL
- Drip tray
- Q-Tight® closing: Special closing to prevent leakage from the top of the bag
- Brushless motor
- Variable speed (1.5 - 3 - 4.5 - 6 strokes/s)
- Variable blending time (1 s - 1 h or ∞ / count-down)
- Voltage and frequency: 100-240 V~ 50-60 Hz / 8A
- Fuses: T8A L 250V 5x20mm
- Max power: 1000 W
- International protection marking: IP21
- 3-year warranty (after recording the warranty card)
- Shock absorbers lifetime warranty
- Window door lifetime warranty
- Availability of spare parts: 10 years
- In compliance with: ISO 7218, ISO 6887-1 and FDA BAM (Bacteriological Analytical Manual)
- Manufactured under CE, RoHS, WEEE
- GLP-compliant use (Good Laboratory Practice)
- Designed and made in France

WEIGHT AND DIMENSIONS

- Dimensions (w x d x h): 51.2 x 52.3 x 46.8 cm, net weight: 49.2 kg
- Box (w x d x h): 61 x 57 x 57 cm, gross weight: 53.9 kg

DELIVERED WITH

- 1 drip tray
- 1 power cable
- 1 user's manual
- 1 free pack of blender bags
- 1 certificate of conformity / warranty card

ACCESSORIES

- **Ref 037 000 JumboStand:** Stainless steel ergonomic table for **JumboMix® 3500**.
Dim. (w x d x h): 58 x 51 x 54 cm. Loading height of the **JumboMix®**: 45 cm