

# Anabac<sup>®</sup> natural

Autoclave deodorant



## Fragrance 100% natural origin

During and after the process, **Anabac<sup>®</sup> Natural** refreshes and perfumes the air in the room where the autoclave is located. **Anabac<sup>®</sup> Natural** is available in 3 fragrances: Coconut, Vanilla et Mimosa.



Coconut



Vanilla



Mimosa

interscience




Based on natural extracts, **Anabac Natural®** efficiently eliminates the odours released by autoclaves and brings a pleasant fragrance to your laboratory.


Biodegradable, this new range is more respectful towards lab technicians and the environment.

**Anabac Natural®** is certified ECOCERT according to the Cosmos-standard.


**Anabac® Natural** is compatible with all sterilization and decontamination autoclaves and does not affect the autoclaving process.



During sample sterilization, bad smell can emerge from the autoclave.



To handle this, place 1 capsule - or 2 if the smell is very strong - in the water at the bottom of the autoclave (for water-bottom autoclave) or in a beaker of water (for steam autoclave).



After a sterilization cycle, the capsule melts, allowing the perfume to be released.



**Coconut**



**Vanilla**



**Mimosa**

<b>Reference</b>	321 100	321 200	321 300
<b>Fragrance</b>	coconut	vanilla	mimosa
<b>Box dimensions</b> (w x d x h)	12 x 10 x 8 cm		
<b>Gross weight</b>	100 g		
<b>Pot of</b>	50 capsules		

**Ask us for samples**

Contact us for **Infos | Videos | Demos | Quote** at [info@interscience.com](mailto:info@interscience.com)

and visit [www.interscience.com](http://www.interscience.com)