

BagMixer® 400 P

400 mL lab blender

Ref 021 230



Standards:



Manufactured under:



Product made for INTERSCIENCE by Interlab, an ISO 9001 certified company

BagMixer® 400 P is an easy-to-use and powerful lab blender. Adapted to all kinds of applications and with a guarantee of optimal bacterial extraction, it is a great tool for lab analyses.

TECHNICAL SPECIFICATIONS

- All stainless steel 304L
- Plain door
- Useful volume: 50 - 400 mL
- Fixed speed (8 strokes/s)
- Variable blending time (30 - 210 s or ∞)
- **Q-Tight® closing: tight closing for no more leaks from the top of the bag**
- Fully opening door (270°)
- Also available in **CoolMix®** version

- Voltage and frequency: 110-127V~ / 220-240V~ 50-60Hz
- Fuses: T 3.15A L 240V / T 1.25A L 5x20 mm
- Max power: 135 W

- Dimensions (w x d x h): 27 x 38.5 x 28 cm, weight: 16.6 kg
- Box (w x d x h): 49.5 x 35 x 36.5 cm, weight: 18.6 kg

- 3-year warranty
- **Shock absorbers lifetime warranty**

- In compliance with: ISO 7218, ISO 6887-1 and FDA BAM (Bacteriological Analytical Manual)
- Manufactured under CE, RoHS, WEEE - Products made for INTERSCIENCE by Interlab, an ISO 9001 certified company
- Designed and made in France

DELIVERED WITH

- 1 power cable
- 1 user's manual
- 1 quick user guide
- 1 free pack of blender bags

OPTIONS / ACCESSORIES

- Ref 026 000 **Driptray:** Stainless steel security drip tray for **BagMixer® 400**