

BagMixer® 400 W

400 mL lab blender

Ref 022 230

WINDOW DOOR
LIFETIME
WARRANTY

SHOCK ABSORBERS
LIFETIME
WARRANTY

3
YEAR
WARRANTY



Standards:



Manufactured under:



Product made for INTERSCIENCE by Interlab, an ISO 9001 certified company

BagMixer® 400 W is an easy-to-use and powerful lab homogenizer. Adapted to all kinds of applications and with a guarantee of optimal bacterial extraction, BagMixer® 400 W, with its window door, ensures permanent control over blending process for perfect results.

TECHNICAL SPECIFICATIONS

- All stainless steel 304L
- Window door
- Useful volume: 50 - 400 mL
- Fixed speed (8 strokes/s)
- Variable blending time (30 - 210 s or ∞)
- **Q-Tight® closing: tight closing for no more leaks from the top of the bag**
- Fully opening door (270°)
- Security drip tray

- Voltage and frequency: 110-127V~ / 220-240V~ 50-60Hz
- Fuses: T 3.15A L 240V / T 1.25A L 240V 5x20 mm
- Max power: 135 W

- Dimensions (w x d x h): 27 x 38.5 x 28 cm, weight: 16.6 kg
- Box (w x d x h): 49.5 x 35 x 36.5 cm, weight: 18.8 kg

- 3-year warranty
- **Shock absorbers lifetime warranty**
- **Window door lifetime warranty**

- In compliance with ISO 7218, ISO 6887-1 and FDA BAM (Bacteriological Analytical Manual)
- Manufactured under CE, RoHS, WEEE - Products made for INTERSCIENCE by Interlab, an ISO 9001 certified company
- Designed and made in France

DELIVERED WITH

- 1 security drip tray
- 1 power cable
- 1 user's manual
- 1 quick user guide
- 1 free pack of blender bags