

BagRack® Slide

Storage rack for bags

Ref 221 440



Manufactured under:



BagRack® Slide enables the clean and methodical storage of 12 bags with samples. Perfect to use with **BagFilter® Pipet & Roll**, which allows an easy pipetting, avoids cross contaminations from the pipet, and allows quick and clean disposal of the bags. The **BagRack® Slide** allows the user to slide the bag from the side for insertion and also access for pipetting. This rack is compatible with pre-enrichment of the sample for the detection of pathogens bacteria such as listeria or salmonella.

TECHNICAL SPECIFICATIONS

- Fits 12 bags
- Compatible with:
 - BagFilter® Pipet & Roll for easy pipetting and quick and clean disposal
 - BagFilter® Roll for quick and clean disposal
 - BagClip® and BagFilter® Pipet for easy pipetting
 - BagClip® and standard 400 blender bags
- All stainless steel
- Autoclavable
- Stackable

- Designed and made in France

BagRack® Slide 400

Ref 221 440

- For 12 bags of 400 mL
- Dimensions (w x d x h): 21 x 35.5 x 22.5 cm
- Box (w x d x h): 23.5 x 37 x 24.5 cm, gross weight: 1.6 kg

OPTIONS/ACCESSORIES

- Ref 221 042 **DripTray®**: Stainless steel drip tray for **BagRack® 400** and **BagRack® Slide 400**
Box (w x d x h): 22 x 42 x 4.5 cm, gross weight: 1 kg