

## BagSeal®

Thermal sealing unit

Ref 261 000



Manufactured under:



Product made for INTERSCIENCE by Interlab, an ISO 9001 certified company

BagSeal® allows a wide and tight sealing of all kind of bags (with/without filter). Perfect closure is guaranteed.

### TECHNICAL SPECIFICATIONS

- All stainless steel
- Large sealing: 245 x 5 mm
- Adjustable heating time
- LED light code at the end of the sealing
- For bags 80/100/400 mL
- Pressure: 22 Newton / cm<sup>2</sup>
  
- Voltage and frequency: 100-127 V~ / 220-240 V~ 50-60 Hz
- Max power: 1200 W
  
- Dimensions (w x d x h): 37 x 27.5 x 32 cm, weight: 10.8 kg
- Box (w x d x h): 49.5 x 39.5 x 49 cm, weight: 13.6 kg
  
- Warranty 1 year
- Product made for INTERSCIENCE by Interlab, an ISO 9001 certified company
- Designed and made in France

### OPTIONS / ACCESSORIES

- Ref 261 500 Maintenance kit for BagSeal®: 2 sealing threads + 2 fuses + 2 Teflon strips