

Beakers (DB50)

Beakers for platers

Ref 415 100



Manufactured under:



DB50 beakers are placed in the easySpiral[®] automatic plater and are containers for the sample (5mL) for dilution or plating.

TECHNICAL SPECIFICATIONS

- Capacity: 5 mL
- Gamma ray treated: Gamma 5 to 12 kGy, with certificate
- Material: Polypropylene
- Autoclavable up to 10 times
- Autoclave the beakers in the tube half-open.
- Average shelf life: 25 years
- Box of 1000 beakers
- Box dimensions: (w x d x h): 12 x 10 x 8 cm, gross weight: 0.15 g