

## CONFORMANCE DECLARATION

We, at interscience, certify that **Scan® software** conforms to **CFR Regulations 21 Part 11**, issued by FDA (Food and Drug Administration of USA) which specifies how to handle electronic records and electronic signatures.


Scan® software has the following features:

- Each user is attributed a login name and password.
- Each user gains access to different software function according to its profile (three user profiles).
- Automatic audit trail on all record changes.
- Electronic signatures (Login name) embedded into data records, which authenticates their author.
- User may save data and results into protected file and generate reports

St Nom-la-Bretèche, on February 28<sup>th</sup> 2017



Emmanuel JALENQUES  
Co-CEO



Jean Marie EPITALON  
Software development  
manager