

dataLink™

Petri dish labeling solution

Ref. 410 150



in compliance with:



Traceability:



Manufactured:



The dataLink® pack is a Petri dish labeling solution.

Thanks to the dataLink® software, users can manage and print the information of their choice on a label, which is then stuck onto a Petri dish.

- After incubation, this data, from an easySpiral® plater, configures the automatic counting of colonies on Petri dishes by a Scan® counter (volume, dilution, time, deposit mode, etc.).
- Imported from an external file, this data will generate a series of labels which will be printed according to a template designed by the operator.

DATALINK™ SUPERVISION AND PRINTING SOFTWARE

- User-friendly interface in English
- Complies with 21 CFR Part 11 regulations requiring the control of systems handling digital data, in order to guarantee their integrity and security.
- User account management (connections, authorizations, action history)
- View audit trail or global activity log
- Monitoring of plating parameters processed by easySpiral® and Scan devices®
- Automatic label printing, triggered by easySpiral® plating
- Import of XLSX, CSV or TXT files with data recognition
- Import of raw data by barcode scanning
- Customized labeling according to a print template
- Share Excel™ files with LIMS or external database

DATALINK™ PRINTER

- Thermal transfer printing
- Chuck diameter: 12.7 mm
- Print resolution: 203 dpi
- Label size: 50 x 10 mm, thickness: 0.08 - 0.19 mm
- Pre-peel labelling

- USB 2.0 port: computer connection
- Supply voltage and frequency: 100-240 V ~ 50-60 Hz
- Power: standby < 10 W / nominal 100 W
- Dimensions (W x D x H): 20.2 x 26.4 x 19.6 cm, net weight: 2.05 kg
- Box (W x D x H): 28.5 x 35 x 25 cm, gross weight: 3.5 kg

- Warranty: 1 year, upon registration of warranty card
- Optimum life of thermal transfer film and one roll of self-adhesive labels: 1 year
- Optimum label storage temperature: between 15°C and 25°C
- Made in Germany
- In compliance with EU regulations and Good Laboratory Practices (GLP)

DELIVERED WITH

- 1 cord with power supply
- 1 USB cable
- 2 user manuals
- 1 thermal transfer film
- 1 chuck
- 2 rolls of 2,500 self-adhesive labels, 1 roll mounted on printer
- 1 USB key with dataLink™ software
- 1 barcode reader

COMPUTER REQUIREMENTS

- Processor: Intel i5, 2.8 GHz, quad Core (or higher)
- Operating system: Windows™ 10 or 11 (or higher)
- RAM: 4 GB
- Equipment: 4 free USB ports (USB hub compatible)
- Screen resolution: 1280 x 1024 pixels (or higher)

OPTIONS / ACCESSORIES

- Ref. 410 120 Printer labels (dataLink™): 2 rolls of 2,500 labels for dataLink™ printer® - Label size: 50 x 10 mm - Quantity: 2 x 2500
- Ref. 410 130 Printer film (dataLink™) : Thermal transfer film for dataLink™ printer