

## JumboMix® 3500 P CC


3500 mL lab blender

Réf. 033 110 (100-127 V~ / 220-240 V~)

Réf. 033 230 (220-240 V~)



Standards:



Manufactured under:



JumboMix® 3500 P CC enables the quick homogenizing of large samples. Avoiding all risks of cross-contamination, it is an easy-to-use and powerful lab homogenizer. To improve blending efficiency and lower preparation time, it is equipped with patented adjustable paddles. Adapted to all kind of applications and with a guarantee of optimal bacterial extraction, it is a great tool for lab analyses.

### TECHNICAL SPECIFICATIONS

- All stainless steel 304 L
- Multifunction digital screen
- Plain door
- Fully opening door (270°)
- Q-Tight® closing: tight closing for no more leaks from the top of the bag
- Click & Clean® removable paddles (patented)
- Side-by-side paddle stop
- Also available in CoolMix® version
  
- Adjustable blending power (patented)
- Optimal blending volume: 200 - 3750 mL
- Variable speed (1.5 - 3 - 4.5 - 6 strokes/s)
- Variable blending time (1 s - 1 h or ∞ / count-down)
  
- Voltage and frequency  
Ref. 033 110 : 100-127 V~ / 220-240 V~ 50-60 Hz  
Ref. 033 230 : 220-240 V~ 50-60Hz
- Fuses  
Ref. 033 110 : T 4 A L 250 V / T 2 A L 250 V 5x20 mm  
Ref. 033 230 : T 2 A L 250 V 5x20 mm
- Max power: 400 W
  
- Protection index : 21
- Dimensions (w x d x h): 52 x 47 x 47 cm, weight: 53.2 kg
- Box (w x d x h): 62 x 58 x 58 cm, weight: 61 kg
- 3-year warranty
- Shock absorbers lifetime warranty
- Availability of spare parts : 10 years

# interscience

- In compliance with: ISO 7218, ISO 6887-1 and FDA BAM (Bacteriological Analytical Manual)
- Manufactured under CE, RoHS, WEEE - Products made for INTERSCIENCE by Interlab, an ISO 9001 certified company
- Designed and made in France

## DELIVERED WITH

- 1 power cable
- 1 user's manual
- 1 quick user guide
- 1 free pack of blender bags

## OPTIONS / ACCESSORIES

- Ref 036 000 DripTray: Stainless steel security drip tray for JumboMix®
- Ref 037 000 JumboStand: Ergonomic table for JumboMix® 3500

### PARIS

Phone: +33 1 34 62 62 61  
[info@interscience.com](mailto:info@interscience.com)

### FRANKFURT

Phone: +49 611 7238 7770  
[sales.germany@interscience.com](mailto:sales.germany@interscience.com)

### BOSTON

Phone: +1 781 937 0007  
[sales.usa@intersciencelab.com](mailto:sales.usa@intersciencelab.com)

### SHANGHAI

电话: +86 178 2123 6642  
[sales.china@interscience.cn](mailto:sales.china@interscience.cn)

### SINGAPORE

Phone: +65 6977 7232  
[sales.asia@interscience.com](mailto:sales.asia@interscience.com)

### TOKYO

Phone: +81 3 6712 9715  
[sales.japan@interscience.com](mailto:sales.japan@interscience.com)