

## JumboMix® 3500 VW

3500 mL lab blender

Ref 032 230

WINDOW DOOR  
LIFETIME  
WARRANTY

SHOCK ABSORBERS  
LIFETIME  
WARRANTY

3  
YEAR  
WARRANTY



Standards:

ISO  
7218

ISO  
6887-1

FDA BAM  
Bacteriological Analytical  
Manual

Fabrication :

RoHS  
2002/95/EC

WEEE  
2002/96/EC



Product made for INTERSCIENCE by Interlab, an ISO 9001 certified company

**JumboMix® 3500 VW** allows quick and easy homogenizing of big samples. To improve blending efficiency and lower preparation time, it is equipped with patented adjustable paddles. JumboMix® 3500 VW, with its wide window door, ensures permanent control over your blending, for an always perfect result.

### TECHNICAL SPECIFICATIONS

- All stainless steel 304L
- **Window door**
- Multifunction digital screen
- Useful volume: 200 - 3500 mL
- Variable speed (1.5 - 3 - 4.5 - 6 strokes/s)
- Variable blending time (1 s - 1 h or ∞)
- **Adjustable blending power (patented)**
- **Side-by-side paddle stop**
- **Q-Tight® closing: tight closing for no more leaks from the top of the bag**
- Fully opening door (270°)
- Security drip tray
  
- Voltage and frequency: 100-127V~ / 220-240V~ 50-60Hz
- Fuses: T 4A L 240V / T 2A L 240V 5x20 mm
- Max power: 400 W
  
- Dimensions (w x d x h): 51.8 x 47 x 47 cm, weight: 50.2 kg
- Box (w x d x h): 62 x 58 x 58 cm, weight: 58.2 kg
  
- 3-year warranty
- **Shock absorbers lifetime warranty**
- **Window door lifetime warranty**
  
- In compliance with: ISO 7218, ISO 6887-1 and FDA BAM (Bacteriological Analytical Manual)
- Manufactured under CE, RoHS, WEEE - Products made for INTERSCIENCE by Interlab, an ISO 9001 certified company
- Designed and made in France

### DELIVERED WITH

- 1 security drip tray
- 1 power cable
- 1 user's manual
- 1 quick user guide
- 1 free pack of blender bags

### OPTIONS / ACCESSORIES

- Ref 037 000 **JumboStand**: Ergonomic table for **JumboMix® 3500**

1 / 1

**JumboMix® 3500 VW\_FT\_0717** Pictures and information are not contractually binding. JumboMix® is trademark. INTERSCIENCE reserves the right to change or improve the specifications of its products without notice. Please visit [www.interscience.com](http://www.interscience.com) for current updates and additional information. RCS 950 356 220 Versailles. INTERSCIENCE SARL F78860

[www.interscience.com](http://www.interscience.com)

interscience

30, ch. Bois Arpents - 78860 St Nom - FRANCE  
T: +33 (0)1 34 62 62 61 - F: +33 (0)1 34 62 43 03  
Email: [info@interscience.com](mailto:info@interscience.com)

interscience USA & CANADA

32 Cummings Park - Woburn, MA 01801 - USA  
T: +1 781 937 0007 - F: +1 781 937 0017  
Email: [sales.usa@interscience.com](mailto:sales.usa@interscience.com)

interscience CHINA

上海市徐汇区吴兴路277号锦都大厦588室-200030  
电话: +86 (0)21-64739390 - +86 189 3097 0733  
邮址: [sales.china@interscience.com](mailto:sales.china@interscience.com)

interscience SOUTH-EAST ASIA

The Centropod, 80 Changi Rd - 05-07 SINGAPORE 419715  
T: +65 6909 0825 - M: +65 8118 5178  
Email: [sales.asia@interscience.com](mailto:sales.asia@interscience.com)