

## JumboMix® 3500 W CC

3500 mL lab blender

Ref. 032 110 (100-127 V~ / 220-240 V~)

Ref. 032 230 (220-240 V~)

WINDOW DOOR  
LIFETIME  
WARRANTY

SHOCK ABSORBERS  
LIFETIME  
WARRANTY

3  
YEAR  
WARRANTY



Standards:

ISO  
7218

ISO  
6887-1

FDA BAM  
Biological Analytical  
Manual

Fabrication :

RoHS  
2002/95/EC

WEEE  
2002/96/EC

CE UK  
CA

MADE IN  
FRANCE

JumboMix® 3500 W CC allows to blend large samples up to 375 g. The adjustable blending power (patented) improves the efficiency of bacterial extraction by adapting to the matrix. JumboMix® 3500 W CC is fitted with removable Click & Clean® paddle system (patented) that simplifies cleaning by providing full access to the blending room. The glass door allows quick controls during the blending process.

### TECHNICAL SPECIFICATIONS

- All stainless steel 304 L
- Multifunction digital screen
- Window door
- Fully opening door (270°)
- Q-Tight® closing: tight closing for no more leaks from the top of the bag
- Click & Clean® removable paddles (patented)
- Side-by-side paddle stop
- Security drip tray
- Adjustable blending power (patented)
- Optimal blending volume: 200 - 3750 mL
- Variable speed (1.5 - 3 - 4.5 - 6 strokes/s)
- Variable blending time (1 s - 1 h or ∞ / count-down)
- Voltage and frequency  
Ref. 032 110: 100-127 V~ / 220-240 V~ 50-60 Hz  
Ref. 032 230: 220-240 V~ 50-60 Hz
- Fuses  
Ref. 032 110: T4A L 250 V / T2A L 250 V 5x20 mm  
Ref. 032 230: T2A L 250 V 5x20 mm
- Max power: 400 W
- Protection index: 21
- Dimensions (w x d x h): 52 x 47 x 47 cm, weight: 50.2 kg
- Box (w x d x h): 62 x 58 x 58 cm, weight: 58.2 kg
- 3-year warranty
- Shock absorbers lifetime warranty
- Window door lifetime warranty
- Availability of spare parts : 10 years

#### PARIS

Phone: +33 1 34 62 62 61  
[info@interscience.com](mailto:info@interscience.com)

#### FRANKFURT

Phone: +49 611 7238 7770  
[sales.germany@interscience.com](mailto:sales.germany@interscience.com)

#### BOSTON

Phone: +1 781 937 0007  
[sales.usa@intersciencelab.com](mailto:sales.usa@intersciencelab.com)

#### SHANGHAI

电话: +86 178 2123 6642  
[sales.china@interscience.cn](mailto:sales.china@interscience.cn)

#### SINGAPORE

Phone: +65 6977 7232  
[sales.asia@interscience.com](mailto:sales.asia@interscience.com)

#### TOKYO

Phone: +81 3 6712 9715  
[sales.japan@interscience.com](mailto:sales.japan@interscience.com)

- In compliance with: ISO 7218, ISO 6887-1 and FDA BAM (Bacteriological Analytical Manual)
- Manufactured under CE, RoHS, WEEE - Products made for INTERSCIENCE by Interlab, an ISO 9001 certified company
- Designed and made in France

**DELIVERED WITH**

- 1 security drip tray
- 1 power cable
- 1 user's manual
- 1 quick user guide
- 1 free pack of blender bags

**OPTIONS / ACCESSORIES**

- Ref 037 000      JumboStand: Ergonomic table for JumboMix® 3500