

JumboMix[®] 3500

Large capacity lab blender



Accuracy and heavy duty blending
for large samples

interscience

Blend up to 3750 mL!

JumboMix® lab blender guarantees optimal bacterial extraction by the action of its 2 paddles. Adapted to any volume between 200 mL and 3750 mL, it guarantees the homogenization in 30 to 60 s of all samples (up to 400 g).

Applications: *E. Coli O157:H7*, trichinella (WarmMix® model), *Salmonella spp*, baby milk powder, chocolate, sterility control tests, Cryptosporidium and Giardia parasites, sample pooling (e.g.: 15 samples of 25 g to reach 375 g)...

Powerful and heavy duty

- Up to 3750 mL capacity
- Ideal for your 375 g samples diluted to 1/10th
- Very silent

Ergonomic

- Side-by-side paddle stop: easy insertion of the bag
- All 304L stainless steel trolley, ergonomic height (in option)
- Digital screen

PATENTED

Easy cleaning

- All stainless steel
- 270° fully opening door
- **Click & Clean®** removable paddles

Security

- The handle shuts and starts the blending
- **Q-Tight®** closing: no leaks upwards
- Security drip tray



Guarantee

- Life-time warranty: Shock absorbers and window
- 3 year warranty (after validation of the card)
- Made in France

PATENTED


Adjustable blending power

- Adapt the blending power to your sample
- Up to 22.6 kg thrust: hard samples are blended




Ideal for your applications

- 375 g for *Salmonella spp*
- 75, 325 or 375 g sample for *E. Coli O157:H7*
- 125 g sample for *Listeria* and *Listeria monocytogenes*
- Sample pooling



SAMPLE POOLING



JumboMix® 3500

3 models available:



JumboMix® 3500 P CC

> Essential features

Ref. 031 230 / 031 110



JumboMix® 3500 W CC

> Window door

Ref. 032 230 / 032 110



JumboMix® 3500 WarmMix® CC

> Thermal regulation up to + 50°C

Ref. 033 230 / 033 110



Technical specifications at the back of the brochure →

JumboLine®


Ideal for your preparation of 375 g samples diluted to 1/10th!



DiluFlow® + JumboKit

Gravimetric dilutor

DiluFlow® Elite - Ref. 503 205 • DiluFlow® Pro - Ref. 502 205 • JumboKit - Ref. 505 006



BagFilter®

Lateral filter bags for lab blenders

2000 mL - Ref. 111 200

3500 mL - Ref. 113 510

BagPage®

Full-page filter bags for lab blenders

2000 mL - Ref. 122 200

3500 mL - Ref. 123 010

BagLight®

Bags for lab blenders

2000 mL - Ref. 132 200

3500 mL - Ref. 133 025



JumboStand

Ergonomic trolley

Ref. 037 000



BagOpen® 3500

Bag opener

Ref. 211 350



BagRack® 3500

Storage rack for bags

2 comp. - Ref. 221 352

4 comp. - Ref. 221 350



BagClip® 3500

Closing clips

Ref. 231 350



BagPipet® & BagTips®

Pipet - Ref. 251 091

Tips - Ref. 252 024

JumboMix_B_0323 Pictures and informations in this brochure are not contractually binding. JumboMix®, JumboMix® Click & Clean®, easySpiral Dile®, Scan®, Q-Tight®, ScanStation®, WarmMix® are trademarks. This document is the translation of the French brochure and is not contractually binding. INTERSCIENCE reserves the right to change or improve the specifications of its products without notice. Please visit www.interscience.fr for current updates and additional information. RCS 950 356 220 Versailles. INTERSCIENCE SARL F78860. Printed by NSI.

JumboMix_B_0323 Pictures and informations in this brochure are not contractually binding. JumboMix®, JumboMix® Click & Clean®, easySpiral Dile®, Scan®, Q-Tight®, ScanStation®, WarmMix® are trademarks. This document is the translation of the French brochure and is not contractually binding. INTERSCIENCE reserves the right to change or improve the specifications of its products without notice. Please visit www.interscience.fr for current updates and additional information. RCS 950 356 220 Versailles. INTERSCIENCE SARL F78860. Printed by NSI.

Products made for
INTERSCIENCE by Interlab,
an ISO 9001 certified company.

ScanStation®
Real-time incubator
& colony counter

Phone: +81 3 6712 9715 - Email: sales.japan@interscience.com

(Z)