

## MiniMix® 100 W CC

100 mL lab blender

Ref 013 230

WINDOW DOOR  
LIFETIME  
WARRANTY

SHOCK ABSORBERS  
LIFETIME  
WARRANTY

3  
YEAR  
WARRANTY



Standards:



Manufactured under:



Product made for INTERSCIENCE by Interlab, an ISO 9001 certified company

MiniMix® 100 W CC is the smallest lab blender in the world. Thanks to its digital screen and window door, control easily the speed and time of blending of your MiniMix® 100 W CC and enjoy a ready-to-analyze sample, with no risk of cross-contamination in less than 60 seconds.

### TECHNICAL SPECIFICATIONS

- All stainless steel 304L
- Window door
- Multifunction digital screen
- Useful volume: 5 - 80 mL
- Variable speed (3 - 6 - 9 - 12 strokes/s)
- Variable blending time (1 s - 1 h or ∞)
- Adjustable blending power (patented)
- Click & Clean® removable paddles (patented)
- Side-by-side paddle stop
- Q-Tight® closing: tight closing for no more leaks from the top of the bag
- Fully opening door (270°)
- Security drip tray
  
- Voltage and frequency: 100-127V~ / 220-240V~ 50-60Hz
- Fuses: T 3.15A L 240V / T 1.25A L 5x20 mm
- Max power: 110 W
  
- Dimensions (w x d x h): 20 x 37.5 x 20 cm, weight: 11.8 kg
- Box (w x d x h): 50 x 30 x 33 cm, weight: 13.4 kg
  
- 3-year warranty
- Shock absorbers lifetime warranty
- Window door lifetime warranty
  
- In compliance with: ISO 7218, ISO 6887-1 and FDA BAM (Bacteriological Analytical Manual)
- Manufactured under CE, RoHS, WEEE - Products made for INTERSCIENCE by Interlab, an ISO 9001 certified company
- Designed and made in France

### DELIVERED WITH

- 1 security drip tray
- 1 power cable
- 1 user's manual
- 1 quick user guide
- 1 free pack of blender bags