

## PetriPile®

### Petri dish storage stack

Ref 241 055 (Ø 55 mm) - Ref 241 065 (Ø 65 mm) - Ref 241 090 (Ø 90 mm)  
Ref 241 150 (Ø 150 mm)



Standards:



Manufactured under:



PetriPile® ensures safe handling, transport and storing of Petri dishes.

### TECHNICAL SPECIFICATIONS

- All stainless steel
- Autoclavable
- Stackable
- Available for Ø 55 mm, Ø 65 mm, Ø 90 mm, Ø 150 mm Petri dishes
- In compliance with ISO 7218
- Product made for INTERSCIENCE by Interlab, an ISO 9001 certified company
- Designed and made in France

#### PetriPile® 55

Ref: 241 055

- For Ø 55 mm Petri dish
- Number of dishes: 36
- Dimensions (w x d x h): 27.2 x 21.5 x 13 cm, net weight: 0,79 kg
- Box (w x d x h): 28.5 x 23 x 15 cm, weight: 1.12 kg

#### PetriPile® 65

Ref: 241 065

- For Ø 65 mm Petri dish
- Number of dishes: 36
- Dimensions (w x d x h): 30.4 x 21.5 x 13 cm, net weight: 0.77 kg
- Box (w x d x h): 31.5 x 23.5 x 14.5 cm, weight: 0.98 kg

#### PetriPile® 90

Ref: 241 090

- For Ø 90 mm Petri dish
- Number of dishes: 36
- Dimensions (w x d x h): 37.5 x 26 x 13 cm, net weight: 1.03 kg
- Box (w x d x h): 38 x 27.5 x 15 cm, weight: 1.4 kg

#### PetriPile® 150

Ref: 241 150

- For Ø 150 mm Petri dish
- Number of dishes: 16
- Dimensions (w x d x h): 37.8 x 35.6 x 13 cm, net weight: 0.88 kg
- Box (w x d x h): 39 x 37 x 14.5 cm, gross weight: 1.26 kg